

Pennsylvania Standards Laboratory

Preparation of Large Volume Provers for Calibration

In order to have your prover calibrated, the Pennsylvania Standards Laboratory (PSL) requires that it be cleaned before arriving at our facility. All surfaces (inside and out), drain lines, pumps, and hoses that may be used during the calibration process must be cleaned prior to calibration, so that no petroleum or other product residue remains that would otherwise spill onto the surface of the testing bay and pollute the sewer system.

The PSL requires the customer to complete the ENTIRE bottom portion of this form for each prover that is being submitted for calibration.

Prover Cleanliness/*Equipment Requirements

| | | | |
|---|--|-------------------------------------|-----------------------------------|
| Prover Manufacturer: | | Prover ID (Serial No): | |
| Prover Volume: | | Minor Graduation Size: | |
| Material: | <input type="checkbox"/> Stainless Steel | <input type="checkbox"/> Mild Steel | <input type="checkbox"/> Aluminum |
| Customer Name: | | | |
| Customer Address: | | | |
| Telephone Number: | | Email Address: | |
| | | | |
| *Prover must have leveling jacks (prover mounted or trailer mounted) present and all jacks must operate freely. | | | |
| *Prover must have functioning spirit levels | | | |
| By signing below, I attest that this prover has been cleaned so that there is no petroleum or other product residue remaining on or in the prover, drain lines, pumps or hoses and also that the *required equipment is adequate. | | | |
| Signature: | | Date: | |
| Printed Name: | | | |
| Present this signed form to the Metrologist on the day of calibration. | | | |